

**Kjeldahl Nitrogen, Total (TKN) 1000mg/L Calibration Standard certified reference material**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | RTC-TKN1000-100ML                 |
| Unit            | 100 mL                            |
| Deliverydetails | No Dangerous Good /not restricted |

**Description**

This certified reference material (CRM) is produced and certified in accordance with ISO Guide 34:2009 and ISO/IEC 17025:2005. All information regarding the use of this CRM can be found on the certificate of analysis.