

**Eukitt® quick-hardening mounting medium,for microscopy (SIGMA-ALDRICH)**

Art. ID                      SAF-03989-500ML  
Unit                        500 mL  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

| Text/Information | Analyte/Parameter                         | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------------|--------------|---------------------|------|--------|--------|
|                  | Eukitt quick-hardening<br>mounting medium | [25608-33-7] |                     |      |        |        |